TANDA

elōs™ professional
hair removal system

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Quick Reference Guide

Important safety information

Please read all warnings and safety information in user manual before using mē.

Note: If you have an active tan do not use mē.

Preparation

We recommend to work in a chilled room where the temperature does not exceed 25°C.

1. Clean your skin and wipe the area with a non-alcoholic wipe.
2. Shave hair in treated area or, if you use the shaver or epilator - trim hair that is longer than 10 mm (you may use the shaver cartridge).

Treatment

1. Press the power button **1** to turn the system on.
2. Select the elōs level by pressing the button near the appropriate level based on the attached table.
3. Press the elōs activation button.
4. Press the epilator/shaver activation button to select desired speed (epilator and shaver are sold separately). Visit www.tanda.com.
5. Ensure you are touching the two grip sensors.
6. Hold the applicator in a position ensuring full contact of the RF delivery bars with the skin. Do not press, just gently touch the skin. System will start emitting elōs pulses.
7. Use the system slowly and continuously gliding on the treatment area without stopping in one place. Try not to repeat treatment over the same spot for at least 10 seconds.

Treatable Skin Types

	Skin Type 1 – Pale white skin (always burns, does not tan) elōs Level: High
	Skin Type 2 – Fair skin (burns easily, rarely tans) elōs Level: High
	Skin Type 3 – Darker white skin (tans after initial sunburn) elōs Level: High
	Skin Type 4 – Light brown skin (burns minimally, tans easily) elōs Level: Medium-to-High
	Skin Type 5 – Brown skin (rarely burns, tans dark & easily) elōs Level: Low-to-Medium
	Skin Type 6 – Dark brown or black skin (never burns, dark tan) elōs Level: Low Device can be used by Skin Type 6 only if the model starts with FG70404

NB: If you feel discomfort, treat at a lower level. If tolerable, increase the elōs setting for higher efficacy.

NB: Do not use on tanned or sun burned skin, as tanning may influence the safety of the treatment. Wait until the sunburn or tan fades prior to performing a test patch on an inconspicuous area. Reexamine the area after 3 days to ensure no side effects have occurred.

NB: Do not tan 24 hours after mē treatment.

NB: During treatment, DO NOT hold the applicator tip still in one spot!

When to use mē

Treatment durations:

For effective results, mē should be used according to the following timings.

1 minute per underarm

2 minutes per bikini area

3 minutes per arm

10 minutes per full leg

1.5 minute per shoulder

5 minutes per abdomen

6 minutes per chest

12 minutes per back

Contact the experts at Tanda to help get you started and set up your treatment routine. Our goal is to help you get the best results possible.

email: info-me@tanda.com
www.tanda.com

When you have finished the treatment

1. Place the applicator in its cradle in the base unit.
2. Turn off the elōs button. The cooling fans will continue operating for 4 seconds. Then turn off the mē device by pressing the power button on the base unit.
3. Gently rinse your skin with lukewarm water and pat dry.
4. Gently apply moisturizing cream or gel aloe vera to your skin (if epilator was used, it is recommended to apply aloe vera).

Table of contents

Quick reference guide	
About the mē hair removal system	7
Box contents	7
Important safety information	8
Recommended elōs levels	11
System description	12
System preparation	14
Using your mē	16
mē treatment procedure	21
Using mē on the face	22
When you are done	26
How to replace your elōs cartridge	27
Troubleshooting	28
Realistic expectations	29
Specifications	32
Labels and symbols	33

About the mē hair removal system

Your new mē system is intended for long-term hair removal and hair reduction. The mē system is easy to use without assistance, and in the convenience of your own home. The system is intended for both women and men of all skin types. The system is designed for the removal of hair from all body parts from the neck and below. A special kit exists for treating dark facial hair. The mē system is powered by elōs™ (Electro-Optical Synergy), clinically proven, gold standard technology that uses low-energy pulsed light together with gentle radio frequency (RF).

NB: Hair reduction results may vary.

Please read the entire user guide and watch the instructional video (available at www.tanda.com) before using the system for the first time. It is recommended to keep the guide in a safe place for future reference and to re-familiarise yourself with this user guide prior to each use of mē.

Box content

Upon opening the mē package you will find the following components:

- Base unit and applicator, including:
elōs cartridge and protective cap for epilator/shaver
- Accessories kit, including:
elōs protective cap
Electrical power cord
Cleaning brush
This User guide (including quick reference guide)
Warranty card
- Additional accessories (sold separately).
Some of these accessories may be included in this package. Actual content is mentioned on the outer package):
Epilator
Shaver
elōs replacement cartridge
Facial kit (use of mē on the face is allowed only after reading the facial treatment instructions in this User Guide and following the instructions).

Important safety information

Please read all warnings and safety information before using the mē system. Always follow basic safety precautions when using any electrical product. The system is designed for domestic use only.

When NOT to use the mē system

You should NOT use the mē system if any of the following items apply to you:

- Tanning: if you are tanned, due to a long skin exposure to the sun, or if you have excessively tanned skin from sun, tanning beds or tanning creams, as tanning may influence the safety of the treatment.
- Wait until the sunburn or tan fades prior to performing a test patch on an inconspicuous area. Reexamine the area after 3 days to ensure no side effects have occurred.
- If you are pregnant or trying to become pregnant (conceive) or nursing.
- If you have a pacemaker or internal defibrillator.
- If you have superficial metal or other implants in the desired treatment area.
- If you have a history of any type of cancer, or pre-malignant moles.
- If you suffer from severe conditions, such as cardiac disorders.
- If you have varicose veins in the treatment area.
- If you have moles, large freckles, darker pigmented areas, scars or birthmarks do not treat over them.
- If you have an impaired immune system due to immunosuppressive diseases such as AIDS and HIV, or use of immunosuppressive medications.
- If you have endocrine disorders such as diabetes. The mē system may be less effective if you suffer from Polycystic Ovaries Syndrome.
- If you have any active condition in the treatment area, such as acne, sores, psoriasis, eczema, burns, herpes simplex, wounds or a rash.
- If you have a history of skin disorders, abnormal wound healing, as well as very dry and fragile skin.
- If you have a history of bleeding coagulopathies, or use of anticoagulants.
- If you are using medications, herbs, food supplements, and vitamins known to induce photosensitivity to light exposure at the wavelength used, such as Isotretinoin (Accutane) and /or topical retinoids within the last six months, Tetracyclines, or St. John's Wort within the last two weeks.
- If you are using medications that may stimulate hair growth, such as Minoxidil.

- If you have had any surgical procedure in the treatment area within the last three months or before complete healing.
- If you have Epilepsy.
- If you have applied any type of anesthetics.
- If you are on gold therapy.
- If the room temperature you are working in is over 25°C (device may not operate).

NB: If you are taking any medication or under medical treatment please consult with your doctor before using mē.

Warnings

- During treatment, DO NOT hold the applicator tip still in one spot! This might cause extensive heat build-up and possible skin damage.
- The system is designed for the removal of human body hair. Any other use is expressly prohibited.
DO NOT use the system for removal of hair on your face or head, unless you have thoroughly read the facial treatment instructions in this user guide.
- DO NOT use the mē system on nipples, areola or your genital parts.
- DO NOT use the mē system over tattoos or permanent makeup.
- NEVER immerse the system or the applicator in water or any other liquid.

- DO NOT rinse any components of the system with any kind of liquid.
- DO NOT use the system where it can fall or be pulled into water. If device is dropped in water, unplug it immediately. DO NOT remove the device from water until unplugged. DO NOT use the device again. Contact the Customer Service Centre.
- If dropped or damaged, DO NOT use the system.
- Always disconnect the power cord from the wall electrical outlet after use.
- DO NOT use the system if any of its components show signs of damage.
- DO NOT open the system's High Voltage Inside enclosure under any circumstances.
- DO NOT operate the device with a cracked or broken glass window on the elōs cartridge. High Voltage Inside.
- To avoid accidents, keep the system away from the hair on your head, and also from clothes, brushes, wires, cords, laces, etc.
- DO NOT use the system on children or allow children to use this device. Keep the system out of the reach of children.
- Adult supervision is recommended for teenagers under the age of 18.
- DO NOT use the system if you have limited physical, sensory or mental capabilities.

Possible side effects

Improper use of the mēsystem could result in possible side effects.

Although these reactions are rare and if they do occur, expected to clear within a short time, it is important to understand that every cosmetic treatment may involve some degree of risk, and that you understand and accept the risks and possible complications that can occur with pulsed light combined with RF hair removal systems.

- **The reaction:** redness or swelling of the skin.

How to treat: may appear either at the time of treatment or after. If there is no improvement within the next 24 hours, consult your doctor.

- **The reaction:** colour change to your skin.

How to treat: stop treatment immediately and consult your doctor.

- **The reaction:** itching of the treated area.

How to treat: cool the area. If the itching continues, consult your doctor.

- **The reaction:** blistering or burning of the skin.

How to treat: cool the area immediately and consult your doctor.

- **The reaction:** unusual pain occurs.

How to treat: stop treatment and consult your doctor.

⚠ This device is not intended for use by individuals (including children) whose physical, sensory or mental capacities are reduced or who have no experience with the device or do not know how to use it, unless they are supervised or have received prior instructions from a person in charge of their safety.

⚠ Children should be supervised to ensure they do not play with the device.

⚠ The devices should not be turned on via an outside timer or a separate remote control system.

⚠ **WARNING:** Keep appliance away from water.

Complete power shutdown will be achieved only when the power cord is unplugged.

Recommended elōs levels

Determine your skin type based on the descriptions below and set your mē system's elōs level accordingly. At the beginning of each treatment area, apply one pulse and evaluate comfort.

If using the recommended elōs level causes discomfort, treat at a lower level.

If comfortable, increase the elōs level for higher efficacy.

elōs cartridge calculator

The amount of time that an elōs cartridge can be used may vary depending on the elōs level setting and your treatment method. You can use our on-line calculator at www.tanda.com.

Treatable Skin Types *

	Skin Type 1 – Pale white skin (always burns, does not tan) elōs Level: High
	Skin Type 2 – Fair skin (burns easily, rarely tans) elōs Level: High
	Skin Type 3 – Darker white skin (tans after initial sunburn) elōs Level: High
	Skin Type 4 – Light brown skin (burns minimally, tans easily) elōs Level: Medium-to-High
	Skin Type 5 – Brown skin (rarely burns, tans dark & easily) elōs Level: Low-to-Medium
	Skin Type 6 – Dark brown or black skin (never burns, dark tan) elōs Level: Low Device can be used by Skin Type 6 only if the model starts with FG70404

NB: If you are actively tanning DO NOT use mē.

*For facial hair removal be sure to refer to the Treatable Facial Skin Types Table in this User Guide.

System description

Your mē system is designed to be easy to use and to maintain.

The system is comprised of the following components:

1. Base unit
2. Control panel
3. Applicator
4. elös cartridge
5. Epilator/shaver protective cap (epilator/shaver will be located here if in use)

Epilator and Shaver are sold separately.
Visit us at www.tanda.com

The treatment is done using the elös cartridge and if you like, you can simultaneously use the epilator or shaver cartridges (Sold separately):

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- Intense pulsed light (IPL) output window.
- Radiofrequency (RF) delivery bars.

2. Epilator cartridge

- Rotating epilating mechanism.

3. Shaver cartridge

- Foil covered shaver combined with a trimmer.

1

2


3


System preparation

The following illustrations display the base unit's control panel and the applicator. The following pages will explain how to prepare the system for use.

Getting acquainted with mē's applicator

Getting acquainted with mē control panel


1. Connect the cable to the system and then to the wall electrical outlet. An orange indicator will appear next to the power button.

2. Before turning on the system, the applicator must be placed in the cradle. Press the power button to turn the system on. All indicators will flash for several seconds, and then the indicator next to the power button will turn green. Now the system is ready for use.

3. Select the elōs level you want to use by pressing the appropriate elōs level button: High, Medium or Low; the active elōs level will illuminate a white indicator light. (Refer to the chart on page 11 of the user guide to set the right elōs level for your skin type).

4. The elōs cartridge is pre-loaded with a set number of pulses for your use; the counter starts at 100% and decreases as you use it. The amount of pulses remaining in the elōs cartridge is indicated by the lights on the panel, as shown. Once the elōs cartridge reaches 5% of its capacity, a red light will start blinking, and at 0%, the red light will remain constant, indicating that the elōs cartridge must be replaced and responsibly disposed of. (Please refer to the WEEE explanation on page 33).

Using your mē

Preparation

We recommend to work in a chilled room where the temperature does not exceed 25°C. Before use, make sure that the RF delivery bars and the light output window on the elōs cartridge are clean (these may be cleaned using a soft tissue or non-alcoholic wet wipes).


Your skin should be clean before starting the treatment. There should be no lotion, make-up, perfume, powder, deodorant or bath/shower oil present on the skin in the area to be treated.

Before treatment, wipe the area with a non-alcoholic wipe.

You can choose one of two options:

1. If you do not use the shaver or epilator accessories, you should shave the required area prior to the elōs treatment.
2. Use with either the epilator or shaver accessory (sold separately) simultaneously during treatment with elōs. In this case, there is no need for shaving prior to treatment.

NB: You can use the epilator or shaver if you have hair up to 10mm.

Epilator sold separately.

How to use

1. After choosing the right elōs level for you, press the elōs activation button on applicator.
2. With the thumb and forefinger of your hand, lightly touch the two metal grip sensors, one on each side of the applicator.
3. Apply the tip of the applicator on your skin so that both RF delivery bars are in complete contact with the skin.
4. Once full contact is achieved, you will notice light flashes, indicating that the system is operating. The elōs activation button will automatically turn itself off if no pulses have been emitted for a while.
5. Move the applicator slowly and continuously without stopping in one place so that the light flashes follow each other and cover the entire area. DO NOT return to the same spot for at least 10 seconds.

DO NOT hold the device in one spot. Facial hair should only be removed after reading the facial treatment instructions in this User Guide. Visit us at www.tanda.com.

NB: It is essential to watch the instructional video (can be watched at www.tanda.com) before using the system for the first time.


NB: It is recommended to use the system in an area with good natural light or in a well-lit environment.

NB: When using the device on difficult to treat areas (eg. underarm) treat the area in 2 directions: left/right and up/down motions, for full coverage. Hard-to-reach areas can also be treated by first using either the shaver or epilator, using the protective cap instead of the elōs cartridge, followed by the application of elōs (while placing the protective cap on the shaver/epilator side). Shaver and epilator attachments sold separately.

NB: The first time you use the system with the epilator cartridge, it may feel slightly uncomfortable, as you are removing hairs from the roots of the follicles. This feeling usually subsides after several treatments, even if you have sensitive skin.

NB: We recommend NOT TO LOOK at the light flashes.

Important Note! The elōs cartridge will not operate unless both RF delivery bars are in contact with the skin, and you are touching the two grip sensors on the applicator.

Important Note! Do not block or cover the applicator's cooling inlet during operation to avoid overheating of the applicator.

HOT hot air will come out of the hot air outlet once the system is in use.

Treating with epilator or shaver cartridge (sold separately)

The mē system is easy to use and easy to care for. Your new mē system has the elōs cartridge and epilator protective cap already inserted into the applicator.

If you choose not to use the epilator/shaver - you must make sure that the epilator protective cap (provided) is inserted, to allow system operation. If you choose to use the epilator or shaver- you should remove the protective cap, and insert the epilator or shaver cartridge instead.


After you have chosen your desired elōs level and activated the elōs cartridge: Press the epilator/shaver cartridge activation button once to run the epilator mechanism at low speed; the number 1 indicator will light up. If desired, press it again to run the epilator at high speed; the number 2 indicator will light up.

NB: The shaver cartridge operates only at Speed level 1.

NB: The elōs cartridge is not interchangeable with the epilator or shaver cartridges.

NB: Epilator and Shaver are sold separately. Visit www.tanda.com

To start treatment please follow instructions on page 21.

Indicator	What does it mean	How to react
Control Panel		
Power light is orange	System in standby mode	Press the power button to switch to "on"
Power light is green	System is on	
"Low" light is white	System elös is set to low level	
"Medium" light is white	System elös is set to medium level	
"High" light is white	System elös is set to high level	
elös cartridge light "100%-20%" is white	Capacity of pulses remaining	
Empty light is blinking red	Only 5% of pulses remaining	Prepare replaceable cartridge
Empty red light is on	elös cartridge is empty	Replace cartridge
Applicator		
elös light is white	elös mode active	
Mechanical hair removal - light number 1 is white	Epilator/ shaver is on speed one	
Mechanical hair removal - light number 2 is white	Epilator is on speed two	

mê treatment procedure

choosing your configuration:

We recommends that you integrate your tried-and-tested hair removing habits with your new mê system:

- If you are used to using a razor, you can shave the treated area before using the elôs cartridge only. Maintain the time of treatment and protocol as illustrated.
- If you are used to using an electric shaver, you can now use the mê system with the elôs and shaver cartridges.*
- If you are used to using an epilator, you can now use the mê system with the elôs and epilator cartridges.*
- If you are used to any combination of shaving or epilating in different areas of your body, use the mê system in the same method.

When treating light hair, it is recommended to use the elôs+epilator combination for higher efficacy.

Treatment durations:

For effective results, mê should be used for at least the amount of time indicated for each body area.

1 minute per underarm

2 minutes per bikini area

3 minutes per arm

10 minutes per full leg

1.5 minute per shoulder

5 minutes per abdomen

6 minutes per chest

12 minutes per back

Using mē on the face

What is facial hair removal?

Removing unwanted hair from the face is very different from removing unwanted hair from other parts of the body. Please be sure to read this section thoroughly before starting to remove facial hair with your mē.

Why is hair removal from your face not like hair removal from other body areas?

- The hair on the face differs from other body areas in its physiology. Its growing cycles are different, and its characteristics, such as thickness (vellus hair), color and depth of follicle are different from body hair.
- The face is more sensitive to hair removal due to sensitivity of nerves in some of its areas.
- The skin on the face and neck is more likely to be exposed daily to the sun's UV radiation and therefore may be more sensitive than other parts of your body.
- Treating near your eyes can cause temporary discomfort, and therefore it is required to use protective eyewear and treat only from your cheeks and down.

When you should not use mē on the face

- If your face is skin type 6 (approved for face skin types 1- 5 only).

- If your facial skin is sunburned or tanned, wait until the sunburn or tan fades prior to performing a test pulse on a hidden area. Evaluate the area after 3 days to ensure no side effects have occurred and only then commence treatment.
- If you have undergone any professional facial treatment (such as chemical peel, laser/light/radio frequency treatment or injected dermal fillers) within the past 3 months.
- If you have undergone any facial surgery (such as face lift, facial implant) within the last three months or before complete healing, please consult with your doctor.
- If you are taking any medications or have a history of any skin disorders, consult with your doctor before using mē on the face or any other body area.
- Please refer to general list of exclusion criteria in the beginning of this User Guide.

Precautions

- Treat below the cheeks line only.

- DO NOT use the epilator or shaver accessories when using mē on the face.
- DO NOT sun tan or use sun beds for at least 2 weeks before a treatment.
- Avoid starting a treatment on your face when you know you will be exposing your face to sunlight in the following 24 hours. If you are going into the sun following a treatment take care to shade the treated area from direct sunlight using a hat or parasol.
- After treatment avoid exposing the treated area to the sun for at least 24 hours.
- Apply sunscreen of at least SPF 30 on the treated area for 2 weeks following the treatment.
- DO NOT treat the eye socket or around the eye area.
- DO NOT use for eyebrow shaping! It is important to avoid accidentally treating the eyebrows as any hair loss is likely to be permanent.

- DO NOT treat over moles or birth marks.
- Protect small areas of sensitive skin such as moles, birthmarks etc by using a white eyeliner. For larger areas of sensitive skin, use a piece of card or thick paper to cover the affected area and prevent exposure to the elōs treatment.
- Protect lips and the lip margin using a white eyeliner.
- DO NOT treat tattoos or permanent makeup.
- DO NOT treat the scalp.
- DO NOT treat inside the nostrils or ears as these areas are delicate.
- Avoid treating areas that contain only vellus hair (very thin, fine hair, also called “peach fuzz”), as this may result in undesired hair re-growth.
- In the very unlikely event that you see pigmented changes on your face please stop the treatment and consult with your doctor.
- Do not use perfumed products or make up, for up to 24 hours after a treatment as this may cause slight irritation.
- If unusual pain occurs during treatment, you should stop treatment immediately.
- Temporary skin redness may occur up to 24 hours after treatment. If the symptoms persist more than 24 hours, and/or you feel significant discomfort

and/or encounter any severe side effects, discontinue treatment and consult with your doctor.

How to use (watch video at www.tanda.com)

1. Before treating your facial hair with the mē please cleanse your skin and remove the hair on the surface of the skin. You may remove the hair on your face with: tweezers, wax, shaver etc.
 2. The treatment should be performed in a broad ambient light environment.
 3. Conceal the delicate areas such as lip edges and facial moles with a white eyeliner.
 4. Turn on your mē device.
 5. Select the LOW elōs level on the base unit and press the elōs button on your applicator.
 6. Use protective eyewear (may be included in your package or in the facial kit you acquired).
 7. Apply a single evaluation pulse on a hidden area on your face or neck and observe treated area. If you observe any unusual reaction, stop the treatment immediately.
 8. Hold the upper part of your mē applicator without touching the two grip sensors and gently position it on the area you want to treat.
 9. To allow an elōs flash on the area you want treat, touch momentarily both grip sensors with the thumb and forefinger of your other hand and release after 1 pulse was emitted.
 10. Without holding the grip sensors, reposition the applicator on the next area to be treated and repeat steps 8 and 9.
 11. Repeat each treated area again 1 or 2 times with at least 10 seconds before you go over the same spot.
 12. If you feel comfortable you may change elōs level from low to med or high.
 13. We recommend treating the face every 1-2 weeks and it is possible to treat the face up to every 2 days!
- * Always start with LOW elōs level. If you feel comfortable, you may increase the elōs level for enhanced efficacy.

grip sensors

Tips

- Please take additional care when treating hair above the upper lip because the central area below the nostrils is particularly sensitive.
- Hair growth under the chin can be denser which can lead to light reddening in the skin following a treatment. This redness is temporary and should completely disappear after 24 hours.
- After the treatment we recommend using an oil free moisturizer for the face in order to feel smoothed.
- For men it is not advisable to treat the beard as it might create uneven results, which may produce changes in facial characteristics that may not be desirable.

For more information regarding facial treatment with mē device, please visit www.tanda.com

NB: Results are based on the completion of protocol and optimization treatments may vary depending on the individual's treatment area, hair type and skin type.

When you have finished the treatment

1. Place the applicator in its cradle on the base unit.
2. Turn off the mē device by pressing the power button on the base unit.
3. Gently rinse your skin with lukewarm water and pat dry.
4. Gently apply moisturizing cream or gel aloe vera to your skin (if epilator was used, it is recommended to apply just aloe vera).

When you are done

- Unplug the system from the electrical wall outlet.
- Unless elōs cartridge is empty, do not remove it from the applicator.
If you do remove the elōs cartridge, insert the elōs protective cap in place.
- If applicable, clean the epilator/shaver cartridge with the cleaning brush (supplied).
It is recommended to clean also with soft tissue or non-alcoholic wet wipes.
- Wipe off the applicator, and the glass window of the elōs cartridge in particular, with soft tissue or non-alcoholic wet wipes.
- Store the applicator back in the base unit.
- Store system away from the reach of children.

Warning: do not rinse any part of the system with water.

NB: Do not use any cleaning agents to clean the mē device and/or its accessories.

Cleaning the shaver cartridge

1. Turn the system off and unplug the electrical power cable from the wall socket.
2. Remove the shaver cartridge by pressing the two release buttons, and sliding it out of the applicator.
3. Remove the foil-covered blade assembly from the shaver cartridge by pulling it out.
4. Use the supplied cleaning brush to clean the blade assembly.
It is recommended to clean also with soft tissue or non-alcoholic wet wipes.
5. Clean the bottom of the shaver cartridge with the cleaning brush supplied.
6. Insert the blade foil assembly into the shaver cartridge again.
7. Put the shaver cartridge back on again. Push down the shaving head until it “clicks” into place.

How to replace your elōs cartridge

1. Make sure that the mē device is turned off.
2. Remove the old elōs cartridge from the applicator by pressing the 2 release buttons on its side and gently pulling it out.
3. Take the new elōs cartridge out of its protective container (avoid touching the glass). Slide the cartridge gently into the applicator until you feel it "click" into place.
4. Before turning on the system, the applicator must be placed in the cradle. Turn on the mē device.
A new elōs cartridge will show 100% on the panel.
5. Turn on the elōs by pressing the elōs button on the side of your applicator. Make sure that the elōs indicator on the applicator is on and is not blinking.
6. Your mē device is ready for use.

Troubleshooting

indicator	What does it mean	How to react
Control Panel		
No lights are on.	Power cable is not connected.	Connect power cable.
Error light is blinking.	<ol style="list-style-type: none"> 1. Button pressed for more than 10 seconds. 2. Grip sensors are held during the start up test. 3. No cartridge level indication. 	<ol style="list-style-type: none"> 1. Release button (system will return to standby mode). 2. Place applicator in cradle. Press power on again. 3. Pull out cartridge and re insert. If not working, insert a new cartridge.
elös level lights are blinking + IPL light pulses stop.	Ambient temperature is too high.	Press elös button again. If not working, move to a cooler area with less than 30 Celsius/86. Fahrenheit degrees.
Continuous error light.	Need to reset the system.	Reset. Disconnect power cord. Wait until all lights are off (up to two min). Reconnect the power cord and turn on the system.
No elös pulses.	Not touching grip sensors. RF bars are not in contact with the treated area. elös button is not pressed.	In order for the device to emit elös pulses, make sure all 3 conditions (elös button on, touching grip sensors and RF bars in contact with the skin) are met.
Applicator		
elös white light is blinking.	elös slot is empty.	Insert elös protective cap or elös cartridge.
epilator/shaver lights are blinking.	Epilator/shaver slot is empty.	Insert epilator, shaver or a protective cap.

If none of these suggested measures seem to help, contact your local provider for support, or contact us at www.tanda.com.

Realistic expectations – what is permanent hair removal with elōs at home?

When an elōs pulse is emitted from the mē device the energy is deposited in the skin, selectively targeting the hair follicle. If the energy reaching the follicle is high enough it can permanently damage the hair follicle and that specific follicle that was damaged will not be able to grow hair any more. Not all hair follicles in the skin area treated by the elōs pulse can be successfully damaged in a given treatment since hair follicles are sensitive to thermal damage only at certain phases of their life cycle. Permanent hair reduction means that when a follicle is thermally damaged it is not able to grow a terminal hair again. Because there are many millions of follicles that can potentially grow hair on the body, and since a significant portion of them is not at the right growth phase, it is not realistic to expect that every single follicle in the treatment area will be disabled. Moreover, many of the follicles in the skin don't have hair that is visible on the surface of the skin. These follicles have a much lower level of sensitivity to thermal damage so that even when they are exposed to the elōs pulse they will not be disabled. These follicles have the potential of growing “new

hair” later in time. Over several treatments with appropriate intervals, most of the follicles that cycle through hair growth will be affected by the elōs procedure and will be affected so they are permanently disabled.

It is important to remember that the body is converting these inactive or dormant follicles to grow later in life just as normal part of hair growth cycle or as a result of many factors: hormonal imbalance, pregnancy, menopause, disease or medications, obesity, or idiopathic causes. This is why there can be no guarantees or promises of “all the hair gone forever” from a certain area. Follicles are closely intertwined with the skin and body response to internal physiological changes and follow a cyclical growth cycle.

Every follicle destroyed will not grow back – your task as a mē user is to treat the majority of follicles with your mē. Follow our treatment protocol and continue treating yourself with mē until you are satisfied with the results.

See that you are treating each area according to the minimum times displayed in the this User Guide and continuously and slowly glide on your skin following the elōs pulses to treat as many follicles possible.

Verify that you glide the applicator in a slow enough motion, so flashes follow each other and are not sporadically emitted on different spots of treated area.

Will I see a permanent result immediately?

Over several treatment sessions there will be less hair growth. This is a gradual and cumulative reduction. When hair does grow back between each of the treatments, there will be less of it, and it will be finer, lighter and slower growing.

With the mē you have the privilege of treating yourself at the privacy of your own home and being able to treat the area continuously when you can see new hair visually appearing in the right growth stage, as a result of the follicle being in the right growth phase. With the mē at home you can follow the recommended protocol and continue optimization treatments as hair appears without having to await an upcoming appointment – but making the appointment with yourself whenever you want and can!

Will ALL of my hair be PERMANENTLY removed from each area treated?

NO. From the reasons mentioned above you cannot expect all the hair to be removed permanently but you will see a significant

reduction in the number of hair that grow on the treated areas. The process is a statistical process and as you progress through the protocol and the optimization treatment you disable a larger percentage of follicles. The results vary from user to user and from one area to another. Because there are many millions of follicles that can potentially grow hair on the body, it is not realistic to expect that every single follicle in the treatment area will be disabled.

This is why there can be no guarantees or promises of “all the hair gone forever” from a certain area. Follicles are closely intertwined with the skin and body, and respond to internal physiological changes.

Treatment tracking

WOMEN

Date: _____

Body Parts:

<input type="checkbox"/>	
<input type="checkbox"/>	
<input type="checkbox"/>	
<input type="checkbox"/>	
<input type="checkbox"/>	

elōs level:

<input type="checkbox"/>	High
<input type="checkbox"/>	Medium
<input type="checkbox"/>	Low

Date: _____

Body Parts:

<input type="checkbox"/>	
<input type="checkbox"/>	
<input type="checkbox"/>	
<input type="checkbox"/>	

elōs level:

<input type="checkbox"/>	High
<input type="checkbox"/>	Medium
<input type="checkbox"/>	Low

Date: _____

Body Parts:

<input type="checkbox"/>	
<input type="checkbox"/>	
<input type="checkbox"/>	
<input type="checkbox"/>	
<input type="checkbox"/>	

elōs level:

<input type="checkbox"/>	High
<input type="checkbox"/>	Medium
<input type="checkbox"/>	Low

Date: _____

Body Parts:

<input type="checkbox"/>	
<input type="checkbox"/>	
<input type="checkbox"/>	
<input type="checkbox"/>	

elōs level:

<input type="checkbox"/>	High
<input type="checkbox"/>	Medium
<input type="checkbox"/>	Low

Date: _____

Body Parts:

<input type="checkbox"/>	
<input type="checkbox"/>	
<input type="checkbox"/>	
<input type="checkbox"/>	
<input type="checkbox"/>	

elōs level:

<input type="checkbox"/>	High
<input type="checkbox"/>	Medium
<input type="checkbox"/>	Low

Date: _____

Body Parts:

<input type="checkbox"/>	
<input type="checkbox"/>	
<input type="checkbox"/>	
<input type="checkbox"/>	

elōs level:

<input type="checkbox"/>	High
<input type="checkbox"/>	Medium
<input type="checkbox"/>	Low

Date: _____

Body Parts:

<input type="checkbox"/>	
<input type="checkbox"/>	
<input type="checkbox"/>	
<input type="checkbox"/>	
<input type="checkbox"/>	

elōs level:

<input type="checkbox"/>	High
<input type="checkbox"/>	Medium
<input type="checkbox"/>	Low

Date: _____

Body Parts:

<input type="checkbox"/>	
<input type="checkbox"/>	
<input type="checkbox"/>	
<input type="checkbox"/>	

elōs level:

<input type="checkbox"/>	High
<input type="checkbox"/>	Medium
<input type="checkbox"/>	Low

Date: _____

Body Parts:

<input type="checkbox"/>	
<input type="checkbox"/>	
<input type="checkbox"/>	
<input type="checkbox"/>	
<input type="checkbox"/>	

elōs level:

<input type="checkbox"/>	High
<input type="checkbox"/>	Medium
<input type="checkbox"/>	Low

Date: _____

Body Parts:

<input type="checkbox"/>	
<input type="checkbox"/>	
<input type="checkbox"/>	
<input type="checkbox"/>	

elōs level:

<input type="checkbox"/>	High
<input type="checkbox"/>	Medium
<input type="checkbox"/>	Low

MEN

Specifications

Light spectrum

- 550 – 1200 nm

Max elos energy 9J/cm²

- Max Optical 4J/cm²
- Conducted RF 5J/cm²
(in equivalent optical units)

Pulsed light repetition rate

- Up to 2 pulses/second (2 Hz)

RF Wave frequency

- 6.78 MHz

Max Radiofrequency power

- 5 Watts

Coverage rate

- Up to 6 cm²/second

Electrical power input

- 100-240 VAC, 50-60 Hz

Dimensions and weight

- 20.5 x 15 x 9 cm [W x L x H]
- 1.17 Kg

Acceptable Environmental Conditions

Storage and transportation

Temperature	-10°C to 55°C
Humidity	90% at 55°C (137°F)
Pressure	700-1060hPa

Operation:

Temperature	10°C to 30°C
Relative Humidity	up to 80% at 37°C
Altitude	up to 3900 m (13,000 ft)

NB: The system will not operate above 30°C.
If device was stored in high humidity, let it dry for an hour before usage

Labels and symbols

The nameplate label is adhered to the bottom of your mē system. It displays the manufacturer's name & address, the system's part & serial numbers, electrical requirements and regulatory standards that the system is designed to comply with.

Do not use this product without first consulting the manual.

Class II equipment

CE mark

WEEE Explanation

This marking indicates that this product should not be disposed with other household wastes throughout the EU. To prevent possible harm to the environment or human health from uncontrolled waste disposal, recycle it responsibly to promote the sustainable reuse of material resources. To return your used device, please use the return and collection systems or contact the retailer where the product was purchased. They can take this product for environmental safe recycling.

Applicable Standards

Safety Standards

- IEC 60601-1
- EN 60335-1: 2002(A1:2004; A2:2009)
- EN 60335-2-8
- EN 60335-2-23: 2003
- EN 60335-2-27: 2003
- EN 60825-9

EMC Standards

- IEC 60601-1-2
- EN 61000-6-1:2007
- EN 61000-6-3:2007

mē is manufactured by Syneron Medical Ltd. Syneron Medical Ltd. reserves the right to make changes to its products or specifications to improve performance, reliability, or manufacturability. Information furnished by Syneron Medical Ltd. is believed to be accurate and reliable. However, Syneron Medical Ltd. assumes no responsibility for its use. No license is granted by its implication or otherwise under any patent or patent rights of Syneron Medical Ltd.

No part of this document may be produced or transmitted in any form or by any means, electronic or mechanical, for any purpose, without the express written permission of Syneron Medical Ltd.

Syneron Medical Ltd. has patents and pending patent applications, trademarks, copyrights, or other intellectual property rights covering subject matter in this document. The furnishing of this document does not give any license to these patents, trademarks, copyrights, or other intellectual property rights except as expressly provided in any written agreement from Syneron Medical Ltd. Specifications are subject to change without notice.

Be sure to read all instructions carefully before starting treatment!

Note: System design may change without prior notification.

Copyright © 2011 Syneron
All rights reserved.

Contact information for customer support:

Syneron Beauty

Website: www.tanda.com,

E-mail: info-me@tanda.com

For Canada only: register your device warranty at tanda.com

Your mē product is warranted to be free from defect in material and workmanship for a period of one year after its original purchase date under normal use. This warranty extends only to the original retail purchaser with original proof of purchase and only when purchased from an authorized Tanda retailer or reseller.

If the product should become defective within the warranty period, contact your local distributor (see www.tanda.com) Syneron Beauty Inc.. In North America, call our customer service team at 519.651.1177 or toll free at 1.877.651.1177 for repair or replacement. If you are outside of North America, contact your local distributor. Contact information can be found at tanda.com. Syneron Beauty reserves the right to replace a defective product with the most comparable product currently available. Your exclusive remedy for any defective product covered by this limited warranty is limited to the repair or replacement of the defective product. Syneron Beauty's entire liability for any defective product shall in no event exceed the purchase price for the defective product. This warranty does not cover products determined by Syneron Beauty to be damaged by any of the following:

- Accident, misuse, abuse or alteration
- Servicing by unauthorized persons
- Use with unauthorized accessories or use other than as instructed
- Connecting to incorrect current and voltage
- Any other conditions beyond our control

TO THE EXTENT PERMITTED BY LAW, SYNERON BEAUTY INC. SHALL NOT BE RESPONSIBLE FOR ANY DIRECT OR INDIRECT DAMAGES RESULTING FROM THE USE OR MISUSE OF THIS PRODUCT INCLUDING ANY GENERAL, INCIDENTAL, SPECIAL, OR CONSEQUENTIAL DAMAGES. ALL IMPLIED REPRESENTATIONS, WARRANTIES OR CONDITIONS INCLUDING BUT NOT LIMITED TO ANY IMPLIED REPRESENTATIONS, WARRANTIES OR CONDITIONS OF FITNESS AND MERCHANTABILITY ARE LIMITED IN DURATION TO ONE YEAR FROM DATE OF ORIGINAL PURCHASE.

This warranty gives you specific legal rights, and you may also have other rights that vary depending on your state, province or other jurisdiction. Some jurisdictions do not allow the exclusion or limitation of incidental or consequential damages, or limitations on how long an implied warranty lasts, so the above limitations and exclusions may not apply to you.

CUSTOMER SERVICE

For questions regarding mē please look for your local service provider at our website www.tanda.com or contact Customer Service.

Customer service contact:

Email: info-me@tanda.com

Hours of operation:

Toll free (NA): 1.877.651.1177

Mon-Fri, 9am-5pm Eastern Standard Time

Syneron Beauty Inc.

6-380 Jamieson Parkway

Cambridge ON, N3C 4N4

THIS HOTLINE CAN BE USED IF YOU
HAVE ANY QUESTIONS REGARDING
USE OF THE DEVICE

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© 2012. All rights reserved. Syneron, Tanda, mē, my elōs, and elōs are trademarks of Syneron Medical Ltd. and may be registered in certain jurisdictions. elōs (electro-optical synergy) is a proprietary technology of Syneron Medical. Patent pending. PB72935